

Device with a daily medicament depot for low-pain injection**Publication number:** DE3604826 (A1)**Also published as:**

DE234613 (A1)

Publication date: 1986-10-16**Inventor(s):** KURZE PETER DR SC NAT [DD]; KRYSMANN WALDEMAR DR RER NAT [DD]; BOCK MANFRED DR MED [DD]; FAEHRMANN KLAUS [DD]; ROEDER HERBERT [DD] +**Applicant(s):** KARL MARX STADT TECH HOCHSCHUL [DD] +**Classification:****- International:** A61M5/315; A61J1/00; A61J1/20; A61M5/178; A61M5/32; A61M5/315; A61J1/00; A61J1/14; A61M5/178; A61M5/32; (IPC1-7); A61M5/20**- European:** A61M5/315D2**Application number:** DE19863604826 19860215**Priority number(s):** DD19850273271 19850214**Abstract of DE 3604826 (A1)**

The invention relates to a device with a daily medicament depot for low-pain injection. The device filled with a medicament can also be worn on the body, has a low weight and is built on the fountain-pen principle. It comprises a dust-proof, thermally insulating casing with an exchangeable, multi-layer, disinfecting, sterile fumigation pad, an inert glass or plastic cylinder with an outside that is firmly gripped, a teflon-coated piston rod and a teflon-coated metal core, as well as a tuned optical and acoustic display system for the injection direction. The device is built on the fountain-pen principle and consists of a syringe with a graduated scale in units of medicament on the unilaterally flattened threaded spindle and a dose syringe with a graduated scale in units of medicament on the unilaterally flattened threaded spindle which coincides with the medicament volume of the syringe. The acoustic system is a bayonet closure with a marking spring extending in a sound space on the unilaterally flattened threaded spindle. The fumigation pad is also used to disinfect the injection site. The puncturing direction can be chosen to be vertical or oblique. The device is particularly suitable for self-injection by diabetics, by diabetics with poor vision or diabetics on holiday or business trips.

Data supplied from the **espacenet** database — Worldwide